


WATER CISTERNS


1 CLEAR ROOF & GUTTERS AND EMPTY CISTERN




2 MAKE DISINFECTING SOLUTION



3 USE DISINFECTING SOLUTION TO SCRUB CISTERN. RINSE WITH PLAIN WATER.




4 FILL CISTERN WITH WATER AND CHLORINE BLEACH. LET SIT FOR 12 HOURS.



5 EMPTY CISTERN

6 REFILL CISTERN



6 RUN UNTIL NO BLEACH SMELL FROM WATER RUNNING OUT. FILL CISTERN.

7 ADD A CHLORINE RESIDUAL



8 MONITOR TO ENSURE A CHLORINE RESIDUAL

